










Overview — Joint Angles





Bilateral Comparison



Source: 22-JUL-2026_pilot_MW_gait_impairment_right
FPS: 59.9
Cycles L: 5 R: 4
Model: mediapipe
Correction: 0.8

Clinical Statistics



Temporal Parameters


Parameter	Value
Cadence	81.8 steps/min
Stride (mean)	1.466 s
Stride (SD)	0.067 s
Step (mean)	0.729 s
Stance L	68.6%
Stance R	60.5%
Double support	29.1%
Cycles	9

ALERTS


* Asymmetry knee: SI=30.6% (>20%)

* Asymmetry ankle: SI=30.6% (>20%)

Trunk and Pelvis Analysis



Pelvis tilt — Not available (lateral view)




Pathology Detection


Pattern	Side	Severity	Value	Confidence	Description
Spastic	Right	mild	33.4	33%	Reduced swing knee flexion (33.4 deg, normal: 60-70)

Normalized Cycles — Left (n=5)


Hip — ROM mean: 27.6°




Knee — ROM mean: 45.0°



Ankle — ROM mean: 21.2°




Trunk — ROM mean: 8.9°




Normalized Cycles — Right (n=4)


Hip — ROM mean: 23.2°




Knee — ROM mean: 33.1°



Ankle — ROM mean: 28.8°




Trunk — ROM mean: 9.3°




Normative Comparison


Hip — vs Normative (adult)




Knee — vs Normative (adult)




Ankle — vs Normative (adult)



Trunk — vs Normative (adult)




GVS / MAP Profile




Quality Dashboard


Detection Rate per Landmark



Per-frame Confidence



Gap Locations (0 gap frames)



Overall Quality Score

96

Quality: GOOD

Detection rate: 100.0%
Mean confidence: 0.88
Gap percentage: 0.0%
Jitter score: 0.0016

Detailed Report

FULL REPORT – GAIT ANALYSIS

Source: /Users/adanmoallemi/Documents/MyoGait/videos/22-JUL-2026_pilot_MW_gait_impairment_right.MP4
FPS: 59.94005994005994 | Resolution: 1920x1080
Model: mediapipe | Correction: 0.8

SPATIO-TEMPORAL PARAMETERS

Cadence: 81.8 steps/min (normal: 100-120)
Stride (mean): 1.466 +/- 0.067 s
Step (mean): 0.729 +/- 0.15 s
Stance L: 68.6% Stance R: 60.5% (normal: ~60%)
Double support: 29.1% (normal: ~20%)

JOINT RANGES OF MOTION (ROM)

LEFT:
Hip -20.6 a 7.0 (ROM: 27.6)
Knee 0.5 a 45.5 (ROM: 45.0)
Ankle -9.4 a 11.7 (ROM: 21.2)
RIGHT:
Hip -12.4 a 10.8 (ROM: 23.2)
Knee 0.4 a 33.4 (ROM: 33.1)
Ankle -12.5 a 16.3 (ROM: 28.8)

SYMMETRY

hip_rom_si 17.3%
knee_rom_si 30.6%
ankle_rom_si 30.6%
step_time_si 1.2%
overall_si 18.5%
(SI < 10% = excellent, 10-20% = good, >20% = asymmetry)

VARIABILITY

Cycle duration CV: 4.8%
Stance % CV: 7.4%

ALERTS

* Asymmetry knee: SI=30.6% (>20%)
* Asymmetry ankle: SI=30.6% (>20%)

Normal reference values:

Hip ROM: 40-50° | Knee ROM: 60-70° | Ankle ROM: 30-40°

Generated by myogait v0.6.0